

SIAction: Non-intrusive, Lightweight Obstructive Sleep Apnea Detection using Infrared Video

You Rim Choi^{*1} Gyeongseon Eo^{*1} Wonhyuck Yoon³ Hyojin Lee¹ Haemin Jang¹ Dongyoon Kim⁴ Hyunwoo Shin^{*2,3} Hyung-Sin Kim^{*1}
Seoul National University¹ Seoul National University College of Medicine² OvarLab³ Columbia University³

INTRODUCTION

- Obstructive sleep apnea (OSA)** is one of the most common sleep disorders




Apnea Hypopnea Index (AHI)

$$= \frac{\text{Apneas} + \text{Hypopneas}}{\text{Total sleep time (hours)}}$$

AHI	Rating
< 5	Normal
5 – 15	Mild OSA
15 – 30	Moderate OSA
> 30	Severe OSA

- Polysomnography (PSG):** Gold standard for sleep evaluation

- ✓ Expensive / Time-consuming / Interscorer discrepancy
- ✓ PSG-based OSA diagnoses can be inaccurate due to
 - (1) First night effect (Attaching sensors in an unfamiliar environment)
 - (2) Single-night stay (No consideration of night-to-night variability)



Before

Appointment & Visit

1. First night & Single night

2. Manual & Time-consuming scoring

3. Diagnosis

After

1. Daily Sleep & Multi-night

2. Auto & Real-time scoring

3. Diagnosis

PRELIMINARY STUDY

Dataset

- Infrared videos from 729 patients, during PSG
- A diverse range of genders and ages spanning from teenagers to individuals over 80 years old
- Various camera angles and room environments across hospitals
- 640 x 480, 5 fps

Advantages

- ✓ Subjects can sleep without sensors; useful for daily sleep monitoring
- ✓ Available anywhere and familiar with users
- ✓ Measuring equipment is not expensive

Challenges




- ✓ Low quality and noisy video
- ✓ Static overall view
- ✓ Not much movement discernible to the naked eyes
- ✓ Blurred information (e.g. subject's face)

- Key Idea using vision AI for OSA detection**

- Difficult to visually distinguish between OSA and normal breathing
- However, in the case of respiratory arousal (RA) accompanied by OSA, the **movement is greater**

Linear relationship between RA frequency and AHI


Indirectly detect OSA by identifying RA !




METHOD

- SIAction** estimates AHI by counting the occurrence of RA during the total sleep time.


Infrared Video Capturing



(1) Input Data Preprocessor




(2) RA Detector




Deep neural network that classifies RA or not
Training data is strategically curated

(3) AHI Estimator



- On-device inference operation**




EVALUATION

Evaluation Dataset

Train	A Hospital, 449 patients			
Valid	A	B	C	Total
	50	0	35	85
	✓ AHI Estimator fitting			
Test	A	B	C	Total
	50	115	80	245

RA Detector

- 9-fold cross validation
- Mean AUC 0.79 (Excellent)



AHI Estimator & OSA Prediction

- High Spearman Correlation coefficient values between the estimated AHI and PSG AHI
- Good OSA prediction performance on hospital A (F1 score 0.876) and B (F1 score 0.891)
- However, the performance for hospital C was relatively lower

→ Refine the AHI estimator using A valid set and C valid set → Shows improved performance

Dataset	Estimator fitting dataset	AHI Estimation Spearman correlation coefficient (ρ)	P-value	Accuracy (%)	OSA Prediction		
A valid (50 patients)	A valid	0.827	1.37e-13	84.0	0.886	0.886	0.886
A test (50 patients)	A valid	0.744	5.89e-10	82.0	0.842	0.914	0.876
B test (115 patients)	A valid	0.756	8.60e-23	83.4	0.867	0.918	0.891
C test (80 patients)	A valid	0.834	8.16e-22	65.0	1.000	0.594	0.745
C test (80 patients)	A & C valid	0.834	8.16e-22	83.7	0.924	0.884	0.903

On-device Performance

Model size (FP16)	Model load Frame capture ready	Frame processing (1 min. clip)	Inference	Total Operation	Peak Memory (RSS)	Peak Memory (Runtime)
5.1 (MB)	1.035 ± 0.007 (s)	0.224 ± 0.042 (s)	3.040 ± 0.046 (s)	3.264 ± 0.088 (s)	839 ± 15.5 (MB)	2.67 ± 0.016 (GB)

CONCLUSION

- We propose a **practical approach** to address PSG limitations for OSA diagnosis, utilizing **human motion** during sleep.
- Our research **pioneers the use of sleep video** to detect OSA events, identifying unique features **with clinical expertise and empirical analysis**.
- We introduce effective **input data curation** for identifying slow and long-term motions in a lightweight manner. **SIAction** is implemented on the **resource constrained device for fast AHI estimation**.
- Through extensive experiments, we demonstrate that **SIAction** maintains **robust accuracy regardless of various environmental changes** such as camera location, sleeping position, blanket occlusion. These results **highlight the broad applicability of SIAction** in diverse real-world settings.

